


# Determination of hydrochloric and hydrofluoric acid in etching baths by thermometric titration



A mix of hydrochloric and hydrofluoric acid can be determined in ethanol- and acetonitrile-containing etching baths by thermometric titration. In the titration curve, two different endpoints can be detected: they are used to quantify the respective acid.

# Method description

## Sample

Simulated etching bath

## Sample preparation

No sample preparation is required

## Configuration

859 Titrotherm	2.859.1010
804 Ti Stand	2.804.0010
800 Dosino, 2x	2.800.0010
20 mL Dosing Unit	6.3032.220
50 mL Dosing Unit	6.3032.250

## Solutions

Titrant	c(NaOH) = 2 mol/L 80 g sodium hydroxide is weighed into a volumetric flask and filled up with deionized water to 1000 mL.
Solvent	1000 mL acetonitrile and 1000 mL ethanol are mixed.

## Analysis

### Blank determination

A linear regression of different sample sizes against consumption is performed. 2.0 mL, 3.0 mL, 4.0 mL, 5.0 mL and 6.0 mL sample solution is pipetted into a titration beaker and 30 mL solvent is added, respectively. The solution is titrated with c(NaOH) = 2 mol/L to one endothermic (hydrochloric acid) and one exothermic (hydrofluoric acid) endpoint.

### Sample determination

The sample analysis is performed in the same way as the blank determination but without the linear regression.

## Parameters

### Blank / Sample determination

Stirring rate	13
Dosing rate	4 mL/min
Filter factor	65

Damping until	0.5 mL
Stop slope	0.200 °C/mL
Stop slope active after	0.5 mL
Evaluation start	0.2
EP criterion 1	55
EP criterion 2	-100
Reaction type 1	endothermic
Reaction type 2	exothermic

## Results

Acid contents (n = 5)

Ratio [HCl:HF]	Recovery HCl / %	S(rel) / %	Recovery HF / %	S(rel) / %
80:20	104.7	0.20	94.8	1.39
60:40	104.3	1.35	104.0	1.69
40:60	104.1	0.87	100.4	2.83
20:80	101.9	1.95	101.5	0.54
10:90	107.2	1.31	101.7	0.21