


Raman-spectroscopic differentiation of inorganic salts



This Application Note shows the ability of handheld Raman systems like the Mira M-1 to identify and differentiate between salts such as carbonates, phosphates, and sulfates. The main focus of this work was to evaluate the influence of the cationic part and the water of crystallization on the Raman-spectroscopic identification of the mentioned salts.

Method description

Introduction

Inorganic salts often have the same anionic part, while the cationic part differs. Many salts can be distinguished by the number of water molecules bond to the salt.

In this study, the influences on the Raman spectrum of the cationic part of the salt and the number of water molecules bound to its anionic part were investigated. Although the differences between the salts are very small, the spectra that were recorded using a handheld Raman spectrometer differed sufficiently to differentiate between them.

Experimental

All spectra were measured using the Mira M-1 Raman spectrometer in auto-acquisition mode, i.e., integration times were determined automatically. A laser wavelength of 785 nm and the Orbital-Raster-Scan (ORS) technique were used. The samples were analyzed in small vials with the vial holder adapter.

Three different anions that give are Raman-active were chosen: carbonate, phosphate, and sulfate. The cationic part of the salts was varied and the resulting spectral changes were investigated.

Results and discussion

Figure 1 shows how the carbonate peak in the Raman spectrum shifts depending on the counter cation. This shift is significant enough to use it to differentiate between the various carbonate salts (except for ammonium carbonate and ammonium bicarbonate, where the difference is very small).

Similar results are also observed for various phosphates; the differences in the spectra are useful for their differentiation (see Fig. 2) – despite the fact that there are some difficulties when trying to distinguish between diammonium phosphate and dipotassium phosphate (due to the similar ionic radius of potassium and ammonium).




Figure 1: Main peaks of the various carbonates.




Figure 2: Main differences in the spectra of the phosphates.




Figure 3: Main differences in the spectra of the sulfates.

Method description

The analyzed sulfates also differed from each other significantly, allowing their unambiguous identification with the handheld Mira M-1, although $\text{CuSO}_4 \cdot 7 \text{H}_2\text{O}$, K_2SO_4 , SnSO_4 , and $\text{ZnSO}_4 \cdot 7 \text{H}_2\text{O}$ could not be identified unambiguously due to their similar spectra.

Conclusions

Generally speaking, handheld Raman spectrometers, such as the Mira M-1, are very useful when it comes to inorganic material identification, even when dealing with salts having the same anionic part. Differences in the cationic parts of the salts help to unambiguously identify many salts with the Mira M-1 handheld Raman analyzer.