

Thermo. Titr. Application Note No. H-089

Title: Automated analysis of hexafluorosilicic acid

Scope: Automated determination of the H_2SiF_6 and HF contents of industrial grade hexafluorosilicic acid

Principle: Titrate a sample of industrial grade hexafluorosilicic acid first with standard NaOH solution to a sharp endpoint, then continue the titration until species have reacted and an excess of NaOH exists. Back-titrate the excess NaOH with standard HCl solution.

The first endpoint of the NaOH titration is used to calculate the "Total" H_2SiF_6 content, which is defined as the "Actual" H_2SiF_6 content plus any residual HF present.

The excess of NaOH is used to force the hydrolysis of SiF_6^{2-} to $\text{SiO}_2(\text{OH})_2^{2-}$ to completion. The back-titration of this excess NaOH with HCl then allows an accurate estimate of the "Actual" H_2SiF_6 content to be made.

Reagents: Titrant 1: 2mol/L NaOH
Titrant 2: 2mol/L HCl
A.R. Potassium hydrogen phthalate

Method: *Basic Experimental Parameters:*

NaOH titration:

Titrant delivery rate (mL/min.)	5
No. of exothermic endpoints	1
Data smoothing factor	75
Maximum dose (mL)	10
Stirring speed (802 stirrer)	8
Delay before start (secs.)	5

HCl titration:

Titrant delivery rate (mL/min.)	5
No. of exothermic endpoints	1
Data smoothing factor	60

Maximum dose (mL)	10
Stirring speed (802 stirrer)	10
Delay before start (secs.)	60
<i>Sample Preparation</i>	
Approximately 1mL of industrial grade hexafluorosilicic acid is weighed into a clean, dry titration vessel. 35mL of D.I. water is added.	
<i>Titration</i>	
A "chained" or "linked" titration program is used to automatically first titrate the acid with NaOH then allow an excess of NaOH to be added before starting a back-titration with HCl to quantify the excess NaOH which has been added. Prior to the commencement of the HCl back-titration, a 60 second delay is programmed to ensure that the hydrolysis reaction has come to completion.	

Example:	<i>Fluka commercial hexafluorosilicic acid, nominal 25%w/w (n=12 replicate titrations)</i>
	"Total" H ₂ SiF ₆ = 24.84±0.07% w/w
	"Actual" H ₂ SiF ₆ = 19.22±0.09% w/w
	Residual HF = 0.78±0.02% w/w
	[F]/[Si] = 6.29±0.01

Calculations:			
$\text{"Total" \% H}_2\text{SiF}_6 = \frac{((V_1 - B_1) \times 144.0918 \times m_1 \times 100)}{(S \times 2 \times 1000)}$			
$\text{"Actual" \% H}_2\text{SiF}_6 = \frac{(V_T - V_1 - (\frac{(V_2 - B_2) \times m_2}{m_1})) \times 144.0918 \times m_1 \times 100}{(S \times 2 \times 1000)}$			
$\% \text{ HF} = \frac{(\% \text{ "Total" H}_2\text{SiF}_6 - \% \text{ "Actual" H}_2\text{SiF}_6) \times 20.00634}{144.0918}$			
$[\text{F}]/[\text{Si}] = \frac{((\frac{\% \text{ H}_2\text{SiF}_6 \times 6}{144.092}) + (\frac{\% \text{ HF}}{20.00634})) \times 144.092}{\% \text{ H}_2\text{SiF}_6}$			
Legend:			
Molecular mass H ₂ SiF ₆	144.0918	Endpoint vol., NaOH titration, mL	V ₁
Molecular mass HF	20.00634	Blank vol. NaOH titration, mL	B ₁
Molarity NaOH	m ₁	Total vol. NaOH titrant delivered, mL	V _T
Molarity HCl	m ₂	Endpoint vol. HCl titration, mL	V ₂
Sample mass, g	S	Blank vol. HCl titration, mL	B ₂

Titration desktop:

