

KF Application Note No. K- 24

Title: Water in smoked fish (salmon, trout)

Summary: The water content of smoked salmon and smoked trout is determined according to Karl Fischer.

Sample: Smoked salmon and smoked trout

Sample

Preparation: Cut up the sample into small pieces. Weigh ca. 2 g sample into a dry flask, add ca. 50 g methanol/chloroform 9 : 1 and mix with the PT 1200 Disintegrator for 120 s. Afterwards stopper the flask with a septum stopper.
For the blank determination a flask is prepared in the same way but without sample.

Instruments and

Accessories: 701 KF Titrino, 720 KFS Titrino or 758 KFD Titrino, 703 Titration Stand, printer, Polytron PT 1200 Disintegrator

Analysis: Pour ca. 25 mL dry methanol into the titration vessel and condition it. Then add ca. 1 g of the prepared sample solution using a syringe and start the water determination.

Reagents:

Solvent: methanol/chloroform 9 : 1,
methanol (dry)

Titrant: Hydranal Composite 5 (Riedel-de Haën)

Results: Salmon: AVG(5) = 63.74 +/- 0.16 % water
Trout: AVG(5) = 66.79 +/- 0.11 % water

Settings: 701 KF Titrino

>titration parameters

extr.time 0 s
stop crit.: drift
stop drift 10 uL/min

>preselections

req.smpl size: on
report: full